


XC-400 Simulator Series



XC-447 Vulva Suturing Training Simulator

Right, left or median posterior incision can be selected to exercise perineum cutting and suturing.

Packing: 8sets/carton, 70x68x52cm, 16kgs

XC-447A Perineum Cutting and Suturing Training Simulator

Feature:

- 1) This model is designed for practice of perineum cutting and suturing
- 2) The main structure of the simulator is a perineum, which is made of simulative foam, and the inner muscle is soft PVC. The simulator is mounted on a plastic framework, easy disassembly.
- 3) There are handles in the back and bottom of the simulator to adjust tension of the perineum.
- 4) Size:36x34x26cm.

Packing: 10pcs/carton, 53x39x47cm, 20kgs



XC-448 Advanced Artificial Abortion Simulated Uterus

- 1) Three pregnant uterus can be opened and add to simulated gestation sac.
- 2) Uterine dilator and curet can be inserted into cervix; can simulate curettage operation and the simulated gestation sac can be scraped.
- 3) With basement support and can fix the uterus in the right position.

Packing: 4pcs/carton, 46x37x27cm, 7kgs

XC-449 Labor Delivery Module

- 1) Six labor stations selected to represent conditions of the cervix and vagina prior to labor, during labor, and at birth
- 2) Stations illustrated are:
STA -5 prior to onset of labor;
STA -4 cervix partially effaced;
STA -3 cervix fully effaced;
STA 0 fetal head at plane of ischial spine;
STA +2 cervix nearing full dilation;
STA +5 crowning of fetal head.

Packing: 1pcs/carton, 43x25x35cm, 7kgs

